
KL 1500 electronic

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GB**Contents****Page**

1. Important hints	8
2. Operation	9
3. Replacing bulb and filters	11
4. Maintenance	11
5. Fault finding	12
6. Accessories for the KL1500 electronic	12
7. Technical data	13

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
1. mains supply
2. on/off switch
3. brightness control
4. light guide collet
5. light guide
6. focusing and filter attachment
7. integrating optics switch
8. bulb holder
9. release catch for bulb housing
10. filter housing
11. halogen reflector bulb
12. air vents
13. air vents
14. safety fuses (on the base of the instrument)
15. model plate (on the base of the instrument)

Cold-light sources are instruments which produce intensive, yet heat-free illumination.

The infrared component of the bulb output is removed by a filter and the remaining visible light is directed via flexible light guides to the object to be illuminated.

1. Important hints

It is important to note that:

- ▶ The working voltage of your KL 1500 electronic must be of the type indicated on the model plate (15).
- ▶ The air vents must be kept free.
- ▶ During operation the light source interior may become hot. The bulb must be allowed to cool down before replacing.
- ▶ The KL 1500 electronic should not be used in a hazardous environment.
- ▶ All repairs must be referred to a qualified electrician.
- ▶ When servicing use only identical replacement parts.

The producer cannot be held liable for any damages caused by non-observance of these instructions.

2. Operation

2.1 Light guide connection

Loosen the light guide collet (4). Insert the light guide up to the stop. For light guides with a locating pin, please ensure that the pin fits into one of the four ferrule slots. Tighten the light guide collet.


2.2 Switching on

On/off switch (2)


Position  : instrument switched off

Position  : instrument switched on

Picture 3: Collet with light guide

Picture 4: On/off switch

2.3 Brightness control

By turning the brightness control (3) the intensity can be adjusted across a wide range. An electronic control ensures a constant luminous output, independent of fluctuations in supply voltage.

The control dial positions are clearly marked enabling repeatable adjustments to be made.

Position  on the scale gives the lowest intensity. Maximum intensity is attained at position 5. The dial is turned past the stop at position  by pushing and continuing to turn. The bulb lifetime is approximately 5 times longer at position  than position 5.

At position 5 the colour temperature is 3,200 K. The dial cannot be turned past position 5.

The introduction of neutral grey filters via the filter holder (10) or the focusing attachment (6) will reduce the intensity without changing the colour temperature.


2.4 Integrating optics

The KL 1500 electronic is fitted with an integrating optics assembly which can be positioned between the bulb and light guide by means of a pivoting lever. This optical assembly creates a more even illumination across the light guide input and will be most effective when the light guide is being used without the focusing lens or filter.


Position  :
Integrating optics positioned in the light path

Position  :
Integrating optics removed from the light path

Picture 5: Brightness control

Picture 6: Integrating optics

3. Replacing bulb and filters

3.1 Replacing the bulb

Open bulb holder (8) by pressing the release catch (9) and pull out to full extent. The current to the bulb is cut off, although the fan will continue to operate. Remove the defective bulb. When inserting the replacement bulb ensure that both levers on the special socket are pushed down. Push the bulb holder back and secure the catch. If the replacement bulb fails to light, switch the instrument off at the on/off switch (2).

3.2 Replacing insert filter

Pull out the filter holder (10) and insert the required filter from the filter set, placing the filter in the holder. By pushing the filter holder to the stop, the filter is positioned in the optical path. The light source can continue to be operated with the filter holder withdrawn.

4. Maintenance


The KL 1500 electronic is virtually maintenance free, with consumable items easily replaced.

Disinfection is not provided for in medical applications. If necessary, the light guides used with the light source have to be disinfected according to the instructions of the producer.

Picture 7: Replacing the bulb

Picture 8: Filter holder

5. Fault finding

Fault	Cause	Remedy
Bulb not lit Fan working	Door not closed	Close door
	Temperature cut-out triggered; Insufficient cooling	Keep air vents clear, avoid high room temperature. Instrument will switch on after approx. 1 min.
	Defective bulb	Replace bulb
Bulb not lit Fan not working	Defective fuse	Replace fuse (14)
	No power	Check power supply
	Defective fuse	Replace fuse (14)

Refer all further repairs to a qualified electrician.

6. Accessories for the KL 1500 electronic

6.1 Light guides

To ensure luminous efficiency use only Schott light guides with your KL 1500 electronic.

An extensive range of accessories is also available; the KL 1500 electronic brochure, which is available on request, gives full details.

6.2 Filters

Optical filters can either be inserted into the filter housing (10) or placed directly at the light guide output by attaching to the focusing lens.

Insert filters

The KL 1500 electronic has a special filter housing which makes it easy to replace even hot filters. The specific filters contained in the insert filter set are detailed in the KL 1500 electronic brochure.

Focusing and filter attachments

The focusing and filter set comprises a focusing attachment with aspherical lens, which fits all light guides in the accessory range (excluding 155 103 with 8 mm dia. fibre bundle), together with 4 screw-fit filters in blue, red, green and orange.

For coloured illumination of larger areas the aspherical lens can be removed. In addition to coloured filters available in the set, standard filters with screw-fit E 19 or push-fit 21 mm dia. may be used.

6.3 Spare parts

Halogen bulb, 15 V/150 W, Typ EFR (Ref.-No. 153 000)

Fuse, for 230 V, primary: (Ref.-No. 153 101)
secondary: (Ref.-No. 153 102)

for 120 V, primary: (Ref.-No. 153 103)
secondary: (Ref.-No. 153 104)

7. Technical Data

Type designation:	KL 1500 electronic
Working voltage:	230 V (Ref. No. 150 100) or 120 V (Ref. No. 150 101)
Frequency:	50...60 Hz (cps)
Power rating:	max. 320 VA
Lamp type:	Halogen reflector lamp 15 V/150 W
Fuse:	for 230 V: primary: T 2 A; secondary: T 10 A for 120 V: primary: T 4 A; secondary: T 10 A
Insulation:	Protection class II Type BF
Dimensions:	Height 143 mm; Width: 183 mm; Depth: 241 mm
Weight:	4.5 kg
Cooling:	Rotary fan
Heat filter:	SCHOTT KG 1.45 x 45; d = 4.2 mm
Emission angle at the end of the light guide:	Without focusing attachment: approx. 70° With focusing attachment: approx. 20°
Max. intensity at light outlet, 2 x dia 4,5 mm: (e. g. 154 202)	Approx. 11 Mlx
Colour temperature at position 5:	Approx. 3 200 K
Colour temperature at position 1:	Approx. 2 800 K
Subject to change	