

EVERIC[®] care

Pharmaceutical glass vial with barrier coating
for high pH fillings with hydrophobic surface



- Established primary packaging solutions reach their limit when protecting sensitive drugs from glass leachables in the high alkaline pH range.
- With EVERIC[®] care, SCHOTT offers a vial with an exceptional ion barrier coating that minimizes leachable levels to the highest extent possible, designed for use in a wide range of pH conditions, even highly alkaline solutions.
- Applied using patented and proven plasma impulse chemical vapor (PICVD) coating technology (known from SCHOTT TopLyo[®] and SCHOTT Type I plus[®] coatings).
- Two additional release criteria have been specifically developed.



Si-O-C-H layer
applied via PICVD



Contact angle for water
> 90° (hydrophobic
surface without silicone)



Long-term
layer stability



Coating bonds covalently
to the glass substrate



Suitable for
depyrogenation



Dense coating
(i. e. non-porous)



Surface shows excellent
barrier properties in
reducing ion leaching,
also in high pH solutions



Stable after
washing process



Compliant with all current
standards, such as Ph. Eur,
USP, JP and CP

Secured drug stability

The excellent barrier reduces ion leaching, providing drug stability even for extreme conditions such as high pH

Improved dosage accuracy
Improved residual emptying and reduced protein adsorption*


Reduced rejects
For lyophilized products, no fogging occurs and an elegant cake is achieved due to the hydrophobic surface

* adsorption behavior needs to be tested case by case

SCHOTT
PHARMA

Release tests

1. Fast performance test



basic

14
13
12
11
10
9
8
7
6
5
4
3
2
1

neutral

acidic

1) Coating attack

Autoclaving filled with 0.005 M KOH (pH 11.7 at start).

2) Coating attack and leaching

Autoclaving filled with 0.1 M HCl (pH 1 at start).

3) Determination of sodium via AAS

Certified release criterion for solution (Na) – limit value defined per format.



2. Drain-off test for hydrophobicity



Autoclaving filled with H₂O.



Half-filled with white milky test fluid and stoppered. Shaken overhead three times.



In specification

Not in specification


Verifications: Superior leachables results for dissolved glass elements in 10R vial with nominal filling volume*

SCHOTT Type I plus® vials served as benchmark. These vials have an SiO₂ inner coating acting as an ion barrier against glass leachables.

Shown here are glass elements in sodium bicarbonate (pH 8.5) after 12 and 48 weeks at 40 °C.

As shown in the figure, for high pH, the barrier properties of EVERIC® care are superior.

* Multiple tests have been performed with different buffer solutions. Please contact SCHOTT for further information



carbon neutral
Schott Pharma
print production

MIX
FSC
FSC C006855

General ordering information

Quality level	TopLine + additional release tests
Packaging	<ul style="list-style-type: none"> • Tray, optionally with divider • Pre-washed & pre-sterilized: adaptiQ®*
Palletizing	Standard Euro pallet (1200 x 800 mm) contains 15–27 layers of nine trays each
Formats	10R*
Pieces per tray	154

* Discuss the availability of other formats with your local SCHOTT representative and the availability of ready-to-use formats (SCHOTT adaptiQ®) – coming soon

Many configurations are available in small quantities as “Fast Track Articles”.

Visit our online shop or speak to your sales representative for more information.

schott-pharma.com/vials

SCHOTT Pharma AG & Co. KGaA, Hattenbergstrasse 10, 55122 Mainz, Germany
pharma@schott.com

SCHOTT
PHARMA