
# KL 1500 HAL User Manual

DE | EN | FR | ES | IT

**SCHOTT**  
glass made of ideas

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
## Instrument overview

- |  |   |
|--|---|
| ① Electronic light intensity setting         | ⑧ Filter slide  |
| ② On/off switch                              | ⑨ Air vents   |
| ③ Locking screw                              | ⑩ Information label                                       |
| ④ Light guide socket with automatic shut off | ⑪ LCD display   |
| ⑤ Carrying handle                            | ⑫ Pivot lever for supplementary optics                    |
| ⑥ Mains connection socket                    | ⑬ Lamp compartment with halogen cold light reflector lamp |
| ⑦ Ventilation grid                           |   |

## Content

1. Important information	20
2. Operation	22
2.1 Light guide connection	22
2.2 Power connection	22
2.3 Start-up procedure	23
2.4 Light intensity setting	23
2.5 LCD display indication and settings	24
2.6 Supplementary optics	25
2.7 Filter slide	26
3. Lamp and socket replacement	28
3.1 Spare Parts	29
4. Maintenance	30
5. Troubleshooting	30
6. Technical data	32

## 1. Important information

### Symbols used

Symbol	Meaning
	Warning of danger (caution, obey documentation)
	Warning of hot surface

### Intended use

The KL 1500 HAL cold light source is intended for use in industrial and laboratory settings<sup>1</sup>. Cold light sources are used to illuminate objects in stereomicroscopy. The visible and intense light from the cold light source is guided to the object via light guides. The light source may only be used with SCHOTT light guides and the SCHOTT accessories intended for this purpose. The light source is only suitable for tabletop use.

### Safety information

Please read and follow these instructions carefully. Ensure that every operator of the device has quick access to these operating instructions.

The protective effect of the device may be impaired if the device is not used as specified by the manufacturer. The manufacturer is not liable for any damage caused by non-observance of these operating instructions.

- Operate the light source only in dry conditions and not in hazardous areas.
- When connecting the device, ensure a safe connection.
- Operate the cold light source only with the voltage specified on the information label.
- Never look directly into the light emission from the light source or the light guide during operation (risk of eye injuries)! The KL 1500 HAL emits visible light of high intensity (according to the EN 62471:2008 standard, risk class 2).
- Never cover the open light guide receptacle or the light guide output (fire hazard)!
- Never cover the open light guide receptacle or the light guide outlet with your hand or other parts of your body (risk of burns)!

<sup>1</sup> You can obtain an individual spectrum of your light source at [info.microscopy@schott.com](mailto:info.microscopy@schott.com)

- When illuminating heat-sensitive or flammable light-absorbing objects, pay particular attention to selecting the light guide spacing and intensity so that no thermal damage occurs to the object. Since light-absorbing materials have the physical property of converting incident light into heat, damage can occur to heat-sensitive or flammable light-absorbing materials.
- All light guide outputs must always be at a safe distance – at least 10 cm – from heat-sensitive or flammable light-absorbing materials when the light source is switched on (avoid possible fire hazard).
- Reduce the brightness and duration of illumination of biological tissue to avoid unnecessary exposure of biological tissue to visible light illumination.
- Keep all ventilation openings free. In case of insufficient cooling, an integrated thermal switch will regulate the light intensity of the light source.
- Safe disconnection from the power supply is only achieved by pulling out the plug.
- Do not open or disassemble the device. Refrain from making technical modifications to the device.
- Allow the filter slider and insert filter to cool down before removing the insert filter. The filter slider must always be in a detent position when operating the light source (see “Filter slide”).
- Allow the light guide to cool down before replacing it. We recommend that you switch off the light source for about five minutes before replacing the light guide.

## 2. Operation

### 2.1 Light guide connection

First, open the locking screw of the light guide holder by turning it counter-clockwise. Then push the light guide in as far as it will go and tighten the locking screw again by turning it clockwise.

When inserting light guides with a guide pin, make sure that it aligns in one of the two guide slots.

As a safety precaution the light guide holder is equipped with a sensor system (automatic switch-off), which allows the lamp to be operated only if a compatible light guide is installed.

If there is no light guide in the light guide socket when the light source is switched on, the error message “No Light Guide” appears in the LCD display.

### 2.2 Power connection

Plug the three-pole cable type H05VV-F3G0.75 (available as an accessory) into the connection socket provided. Only operate the KL 1500 HAL cold light source within the specified voltage range. When replacing the power cord, only use power cords of the same type.

The plug serves as an all-pole disconnect device and must be easily accessible at all times.

The supply cable must not be longer than 3 m.

You can switch KL 1500 HAL on and off by pressing the main power switch. If the device is switched on, the indicator in the LCD display lights up. To disconnect the power supply, remove the power plug. To protect the halogen lamp, the KL 1500 HAL is equipped with a soft start that reduces excessive inrush current.

Adjust the light intensity by turning the light intensity control. The light intensity is shown in the LCD display in values from 0% to 100%. Alternatively, you can also display the light intensity in the equivalent color temperature [Kelvin]. To do this, press the light intensity control once. Please note that this is not an exactly measured color temperature, it only serves as a reference value.

Depending on the light intensity, three different values are shown on the display, which give an indication of the service life of the lamp:

Display	Approximate lamp life
ECO	> 1500 h
STANDARD	150 – 1500 h
BOOST	50 – 150 h

**Attention:** The device is not switched off in the 0 control position, therefore the display lights up. To switch off the KL 1500 HAL, press the power switch.

## 2.3 Start-up procedure

## 2.4 Light intensity setting

## 2.5 LCD display indication and settings

### Display brightness

The brightness set on the light intensity control can be displayed either in percent or Kelvin (see “Light intensity setting”).

Briefly press the light intensity control to switch between the settings.


### Lamp operating time display

The total operating time of the halogen lamp used can be shown in the LCD display in full hours.

Press the light intensity control for longer than five seconds to enter menu control. As soon as the light source is in menu control, the lamp switches off.

Turn the light intensity control to the right and left to move the display arrow in the corresponding menu fields. The respective active menu field is highlighted. Briefly press the light intensity control to enter the selected submenu.

Select the “On/Off” submenu to display the total operating time. Turn the light intensity control to select “On” and confirm by briefly pressing it (check mark appears in the menu field). Select “Menu” and then “Exit”. The lamp switches on again. Select the “On/Off” submenu again to switch off the total operating time. By turning the light intensity control, now select “Off” and confirm it by briefly pressing it (check mark appears in the menu field).

Select “Menu” and then “Exit”. You can reset the total operating time after a lamp change. To do this, select the “Reset” item in the menu control and confirm it by pressing briefly (a check mark appears in the menu field). Select “Menu” and then “Exit”.

Swiveling the additional optics ensures a uniform, high-intensity illumination is achieved even when using light guides with a bundle diameter smaller than 9 mm. If the illumination includes a focusing lens, a uniform illumination can be achieved by swinging out the additional optics. Turn the swivel lever to the left or right to the corresponding symbols to swivel the additional optics in or out. Make sure that the additional optics are always swiveled in or out to the end stop.

## Supplementary optics

## 2.7 Filter slide


The KL 1500 HAL has a filter slide that can be equipped with an insert filter. Various filters are available as accessories.

When inserting the filter slide into the cold light source, there are two detent positions which are signaled by a clicking sound. The filter slide must be in one of the two end positions or the detent position when operating the light source. If the filter slide is in an intermediate position, it may damage the light source.


### Insert filter

Make sure that the filter slide has cooled down before inserting a filter. Wear lint-free gloves to load the filter slide.


1. Pull out the filter slide completely and turn it so that the labeling is legible.

2. Insert the insert filter from the top right behind the two retaining lugs marked in black.


3. Click the insert filter firmly into place by applying slight pressure at the top right. The filter now lies flat in the holder.

To remove, bend the retaining clip slightly to the left and remove the filter.

#### Remove filter

If you want to operate the light source without a filter for a short time, pull the slider out only to the first detent point. In this position, the filter is still in the light source, but no longer in the beam path.

### 3. Lamp and socket replacement

Safety information 

Before changing, switch off the light source and disconnect the power plug. Do not reach into the device immediately after switching it off, as there is a risk of injury from the fan running after it.

Make sure that the lamp and lamp socket have cooled down before replacing them. 

Open the lamp compartment by pressing down on the recess provided and simultaneously pulling it as far as it will go.

#### Lamp replacement

Remove the lamp, including the base, upwards from the holder. Carefully remove the lamp from the base and attach a new lamp (please use the lamp as specified in item 3.1). Then push the lamp into the holder again from above as far as it will go.

Push the lamp compartment in until it locks into place (you will hear it click into place). Now you can plug in the power plug and switch on the light source.

The life of the lamp socket is approximately 20,000 hours of operation or alternatively approximately 20 lamp replacements.

To replace the lamp socket, first carefully pull off the lamp as described above. Use a screwdriver to loosen the upper fastening screw of the clamp and remove the complete component. Then insert the two cable ends of the new base (see “3.1 Spare parts”) into the corresponding openings of the terminal and retighten the screw. Make sure that it is firmly seated.

Put the lamp back on and press it into the holder as described above. Then push the lamp compartment back in until it locks into place (you will hear it click into place). Now you can plug in the power plug and switch on the light source.

Halogen lamp 15 V / 150 W  
Philips 6423  
Osram 64634 HLX

Order No. 153000

Lamp base

Order No. 150710

## Socket replacement

### 3.1 Spare parts

SCHOTT only guarantees proper function, safety and optimum light output when the specified types of spare parts are used.

## 4. Maintenance

Your KL 1500 HAL is maintenance-free.

Use a soft dry cloth or commercially available plastic cleaning cloths for external cleaning of the device.

## 5. Troubleshooting

If you are unable to rectify the problem using the measures described below, please contact your specialist dealer or the nearest SCHOTT agency. After repair, check that the device is safe to use.

Type of fault	Display indication	Possible cause	Remedy
1. Lamp off, fan not running	No display	<ul style="list-style-type: none"> <li>Instrument not switched on</li> <li>Plug not in socket</li> <li>No mains voltage</li> </ul>	<ul style="list-style-type: none"> <li>Switch the instrument on</li> <li>Insert the plug</li> <li>Check mains voltage</li> </ul>
2. Lamp off, fan running	Open Circuit	<ul style="list-style-type: none"> <li>Lamp compartment not connected</li> <li>Lamp defective</li> <li>No lamp</li> <li>No lamp contact</li> </ul>	<ul style="list-style-type: none"> <li>Close lamp compartment</li> <li>Replace lamp (see point 3)</li> <li>Insert lamp</li> <li>Press lamp into the socket, replace the lamp base if necessary</li> </ul>

Type of fault	Display indication	Possible cause	Remedy
3. Lamp off, fan running	Over-heating	<ul style="list-style-type: none"> <li>• Electronics overheated</li> </ul>	<ul style="list-style-type: none"> <li>• Provide for sufficient cooling, verify correct lamp type, place the instrument back into operation after an extended cooling period</li> </ul>
4. Lamp off, fan not running	No Fan	<ul style="list-style-type: none"> <li>• Fan defective, impeller jammed</li> </ul>	<ul style="list-style-type: none"> <li>• Have authorised customer service perform repairs</li> </ul>
5. Lamp off, fan running	No Light Guide	<ul style="list-style-type: none"> <li>• No recommended light guide in the light guide socket</li> <li>• Light guide is not plugged in until it stops</li> </ul>	<ul style="list-style-type: none"> <li>• Plug recommended light guide into the light guide socket, check for proper seating</li> <li>• Plug in light guide into the light guide socket until it stops</li> </ul>


## 6. Technical data

Properties		Values
<b>General information</b>		
Type description		KL 1500 HAL
Dimensions (W x D x H)	mm	approx. 202 x 286 x 146
approx.	kg	approx. 4.2
Cooling		Axial fan
Operating environment*	°C	+ 5 ... + 40 Indoor use only
Relative air humidity*	%	up to 31°C ambient temperature: 85% 31°C to 40°C ambient temperature: decreasing linearly to 75%.
Air pressure*	hPa	800 ... 1060
Max. altitude		≤ 2000 m asl
Transport and storage		
• Temperature	°C	-20 ... +70
• Rel. air humidity	%	10 ... 95 (non-condensing)
• Air pressure	hPa	500 ... 1200
Contamination level		2
IP Protection type		IP 20 (Following IEC 60529)
<b>Electrical information</b>		
Operating voltage, frequency		100 – 240 V ~   50/60 Hz   ±10 %
Power consumption	VA	max. 180
Protection class light source		I
Overvoltage category		II
Lamp type		Halogen reflector lamp
Lamp rated voltage	V	15
Lamp rated power	W	150
Average lamp lifetime		
Level 50%	h	1500
Level 80%	h	150
Level 100%	h	50

\* Test conditions of DIN EN 61010-1 and UL61010-1 standards

Properties	Values	
<b>Lighting information</b>		
Active area of the light guide	mm	9
Total luminous flux at light guide exit (SCHOTT light guide Ø 9 mm, typical value)		
Level 50%	lm	250
Level 80%	lm	450
Level 100%	lm	600
Heat protection filter	mm	SCHOTT KG 2, 45 x 45 d = 2.0 mm, toughened
Light entry angle ( $2\alpha_{\text{eff}}$ )	°	approx. 70
Conformity	CE, UKCA	
EMC emission class	B	

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**WEEE declaration**

Your SCHOTT product was produced and developed with high quality materials and components. The symbol indicates that electrical and electronic devices must be separated from domestic waste and appropriately disposed of after useful life. SCHOTT AG Lighting and Imaging has arranged a waste management system for recycling. Please use this system for removal and help to protect the environment we live in.

Further information regarding our waste management system, please refer to [www.schott.com/lightingimaging/recycle](http://www.schott.com/lightingimaging/recycle)

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SCHOTT AG, Hattenbergstrasse 10,