
# SCHOTT KL 1500 KL 1500-T

SCHOTT

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
1

**SCHOTT Cold light source  
KL 1500/KL 1500-T  
Operating instructions**

- ① mains supply
- ② model plate (on the base of the instrument)
- ③ on/off switch
- ④ air vents (also at back of instrument)
- ⑤ brightness control
- ⑥ halide reflector bulb
- ⑦ door with catch
- ⑧ heat filter
- ⑨ filter holder
- ⑩ light guide socket
- ⑪ focusing and filter attachment
- ⑫ light guide
- ⑬ fuse

**1. Important hints:**


It is important to note that:

- the working voltage of your cold light source KL 1500/KL 1500-T must be of the type indicated on the model plate ②
- the air vents ④ must be kept free
- only SCHOTT light guides guarantee correct function and optimised light output.
- This instrument should not be used in hazardous environment.
- All repairs must be referred to a qualified electrician.
- Bulbs and filters are hot during operation. Allow sufficient cooling time when replacing these parts in order to avoid burnings.

## 2. Operation

### 2.1 Light guide connection

Slacken light guide bush ⑩; insert light guide as far as possible; for light guides with locating pin, please ensure that the pin fits into one of the four ferrule slots; tighten light guide bush.

### 2.2 Switching on/economy setting


On/off switch ③:

Position 0 = instrument switched off  
Position ④ = bulb operates at reduced voltage, which gives a greatly extended operational life.

Colour temperature:  
KL 1500 approx. 3000 K  
KL 1500-T approx. 2900 K

Position ⑤ = operation of bulb at full capacity.

Colour temperature:  
KL 1500 approx. 3200 K  
KL 1500-T approx. 3000 K

### 2.3 Brightness control

Your KL 1500/1500-T offers brightness regulation in graduated stages by means of the brightness control ⑤, without any change in colour temperature. Position 1 of brightness control gives the lowest brightness. Turning the control in a clockwise direction doubles the brightness at each stage. There is a distinct click-stop at each stage. Between stages, the light guide is completely darkened.

Application of neutral grey filters in the filter holder ⑨ or on the focusing attachment ⑩ will extend the lower limit of the control range.

The knob cannot be turned past position 5.

### 3.1 Replacing the bulb

Open door ⑦ by depressing catch. Allow filter holder ⑨ and bulb ⑥ to cool. Release filter holder by turning lever ④, swing filter holder forward. Remove defective bulb from the front. Insert new bulb\*, replace filter holder, ensuring that the pin is properly engaged. Close door.



When the door is open, the current to the bulb is cut off, with continued operation of fan.

**\*Note:** Please check pins, even of new bulbs, for corrosion. In any case, it is advisable to lightly sand the pins to improve the contact.

### 3.2 Replacing heat filter

Open door ⑦ by depressing catch. Allow filter holder to cool. Swing stop lever ⑧ to the left. Remove defective filter and insert replacement into the rear fitting, without touching polished faces with fingers. Remove stains and finger prints immediately with alcohol. Close door.

**Warning: Never operate the KL 1500/KL 1500-T without the heat filter!**

#### 4. Fault finding

Fault	Cause	Remedy
Bulb not lit Fan working	Door not closed	Close door
	Temperature cut-out triggered; Insufficient cooling	Keep air vents ④ free, avoid high room temperature. Instrument will switch on after approx. 1 min.
	Corroded bulb pins	Sandpaper bulb pins to improve contact
	Defective bulb	Replace bulb
Bulb not lit Fan not working	No power	Check power supply
	Defective fuse	Replace fuse ⑬

Refer all further repairs to qualified electrician.

#### 5. Maintenance

The KL 1500/KL 1500-T is maintenance free.

#### 6. Accessories for the

##### KL 1500/KL 1500-T

An extensive range of accessories is available for your KL 1500/KL 1500-T

##### 6.1 Light guides

Designation	Dia. of fibre bundle in each arm (mm)	length (mm)	Ref. No.
Swan-Neck Light Guide 1-arm	4.5	750	154 101
Swan-Neck Light Guide 2-arm	4.5	750	154 201
Swan-Neck Light Guide 3-arm	4.5	750	154 301
Flexible Light Guide	3	1000	155 101
Flexible Light Guide	5	1600	155 102
Flexible Light Guide	8	1000	155 103
Flexible Light Guide 2-arm	4.5	1000	155 202
4-Point Ring Light, inner dia. = 30 mm	4	1000	157 401
4-Point Ring Light, inner dia. = 66 mm	4	1000	157 402
4-Point Ring Light, with magnifying lens, inner dia. = 66 mm	4	1000	157 403
Annular Ring Light inner dia. = 58 mm	9	1000	157 410

## 6.2 Filters

Optical filters are available for fitting to the filter holder ⑨ or with a screw or push fitting for use in combination with the focusing attachment ⑩ at the end of the light guide.

### Insert Filters

Designation	Type	Dimensions (mm)	Ref. No.
Heat Filter	KG 1	42 × 42 × 3	158 300
Set of insert filters comprising 1 each of the following			
Conversion filter	511	40 × 40 × 2	Set: 158100
Blue filter	BG 37	40 × 40 × 3	
Green filter	VG 9	40 × 40 × 3	
Yellow filter	OG 515	40 × 40 × 3	
Red filter	RG 610	40 × 40 × 3	
Neutral Grey filter 1:10	NG 4	40 × 40 × 2	

Alternatively, other hardened filters of optical glass may be used, provided that they measure 40 × 40 and are not thicker than 3 mm.

**On no account may plastic filters be used. The heat filter must not be removed, even when the colour filters are used.**

### Focusing and filter attachments

Focusing and filter set comprising focusing attachment with aspherical lens to fit all swan-neck and flexible light guides (except Ref.No.155103 with 8 mm dia. fibre bundle), and 4 screw-fit filters in blue, red, green, orange

Ref.No. 158 200

Focusing attachment without filters

Ref.No. 158 210

### Polarising filter

with screw-fit for focusing attachment

Ref.No. 158 205

For coloured illumination of larger areas, the aspherical lens of the focusing attachment may be removed.

Apart from the coloured filters available in the set, standard filters with screw-fit M 19 × 0.5 or push-fit 21 dia. may be used.

### 6.3 Additional Accessories

Spare Halogen Bulb 15 V, 150 W

Ref.No. 153 000

Spare Fuse

For 230 V Ref.No. 153 100

For 115 V Ref.No. 153 105

## 7. Technical Data

Type Designation:	KL 1500/KL 1500-T
Working Voltage:	230 V (Ref. No. 150 000) or 115 V (Ref. No. 150 005) 230 V (Ref. No. 150 010) *)
Frequency:	50/60 Hz
Power Rating:	Max. 200 W
Type of Bulb:	Ellipsoid Halogen Reflector Bulb 15 V/150 W
Fuse:	For 230 V: 1.25 A For 115 V: 2.5 A
Insulation:	Protection Class II
Dimensions:	Height: 160 mm Width: 205 mm Depth: 170 mm
Weight:	3.5 kg
Cooling:	Rotary Fan
Heat Filter:	SCHOTT KG 1 (3 mm thick)
Angle of emission at the end of the light guide:	Without focusing attachment: approx. 70° With focusing attachment: approx. 20°
Max. illumination at light guide inlet, dia. 4.5 mm:	Approx. 20 Mlx Approx. 16 Mlx *)
Max. illumination at light guide outlet, dia. 4.5 mm:	Approx. 10 Mlx Approx. 8 Mlx *)
Colour temperature at position $\varnothing$ :	Approx. 3 200 K Approx. 3 000 K *)
Colour temperature at position $\varnothing$ :	Approx. 3 000 K Approx. 2 900 K *)

Subject to alteration without notice

\*) only for KL 1500-T