




# EVERIC<sup>®</sup> amber

Reliable light protection and regulatory compliance



- While every protein is sensitive to light, the more complex the molecule, the higher the possibility of light-induced aggregation or degradation.
- Antibody-drug conjugates (ADCs) are particularly sensitive to light due to their fragility. Even light exposure during filling prior to lyophilization can cause issues.
- Adapting the lighting set-up during filling is cumbersome and may not be sufficient.
- Existing solutions do not comply with USP, Ph. Eur., and JP light transmission requirements, so multiple primary packaging articles may be necessary.
- EVERIC<sup>®</sup> amber offers a solution that provides reliable light protection while remaining compliant with USP, Ph. Eur., and JP standards.



FIOLAX<sup>®</sup> amber glass tubing –  
absorption of light



Improved forming process: tight  
control of dimensional parameters  
and temperature ranges



Unchanged glass composition –  
Type I Borosilicate Glass



Full compliance with light  
transmission requirements of  
Ph. Eur., USP and JP



## Light protection

Protection against ultraviolet rays  
and short-wave visible light



## Global compliance

Suitable for pharmaceutical companies  
operating in multiple regulatory environments

**SCHOTT**  
PHARMA

Release test


01 Mechanical cut & cleaning

- Vial is:
- 1) Cut in vertical direction
  - 2) Cleaned to get a quasi-plan, contaminant-free surface




02 Spectral screening

Light transmission is measured through the middle of the vial body



03 Spectral evaluation

Certified release criteria acc. to wall thickness and filling volume



Verifications: The right production set-up is key to achieving global compliance

Annealing after hot-forming is necessary to avoid residual stress, but it influences transmission.

Light transmission compliance with Ph. Eur., USP and JP standards can be achieved for amber vials with a 1 mm thick wall (2–10R) via tight process control.




Figure 1: Light transmission compliance for 10R vial (wall thickness 1 mm)

By reducing the wall thickness of a 20R vial to 1 mm, compliance with USP, Ph. Eur., and JP standards can be achieved.

Strength-optimized parameters with tight dimensional tolerances, dedicated inspection and packaging with dividers ensure strength similar to conventional vials.




Figure 2: Strength tests for 20R vial with a reduced wall thickness (1 mm)

General ordering information				
Quality level	StandardLine & release test			
Packaging	<ul style="list-style-type: none"><li>Tray (with divider for 10R and 20 ml format)</li><li>Pre-washed &amp; pre-sterilized: adaptiQ® (tray, cup, nest)</li></ul>			
Possible combinations	EVERIC® amber can be combined with EVERIC® Iyo			
Palletizing	A standard Euro Pallet (1200 x 800 mm) contains 15 – 27 layers of 9 trays each			
Formats	2R	6R	10R	20 ml
Pieces per tray	344	186	104	77

Order now:  
[www.schott-pharma.com/vials](http://www.schott-pharma.com/vials)

[schott-pharma.com/vials](http://schott-pharma.com/vials)

SCHOTT Pharma AG & Co. KGaA, Hattenbergstrasse 10, 55122 Mainz, Germany  
[pharma@schott.com](mailto:pharma@schott.com)

SCHOTT  
PHARMA