

The image shows a close-up, low-angle view of several parallel glass tubes in a photobioreactor system. The tubes are supported by blue, ribbed metal brackets. The liquid inside the tubes is a vibrant green, indicating the presence of algae or other photosynthetic organisms. The lighting is bright, creating a high-contrast scene with many reflections on the glass surfaces.

SCHOTT
glass made of ideas

Tubular Glass Photobioreactors

Ensuring efficiency

Pioneering – responsibly – together. These attributes characterize SCHOTT as a manufacturer of high-tech materials based on specialty glass. Always opening up new markets and applications with a pioneering spirit and passion – this is what has driven the **#glasslovers** at SCHOTT for more than 130 years. Represented in 34 countries, the company is a highly skilled partner for high-tech industries: Healthcare, Home Appliances & Living, Consumer Electronics, Semiconductors & Datacom, Optics, Industry & Energy, Automotive, Astronomy & Aerospace. As a foundation company, SCHOTT has anchored responsibility for employees, society and the environment deeply in its DNA. The goal is to become a climate-neutral company by 2030.

With a production capacity of more than 190,000 tons and production sites in Europe, South America and Asia, SCHOTT Tubing is one of the world's leading manufacturers of glass tubes, rods and profiles. More than 60 different glass types are produced in a large variety of dimensional and cosmetic specifications based on a standardized production process and a global quality assurance system. SCHOTT Tubing provides customized products and services for international growth markets such as pharmaceuticals and electronics as well as industrial and environmental engineering.




Contents

| | | | |
|----|--------------------------------------|----|--|
| 4 | ALGAE PRODUCTION SYSTEMS | 16 | BOROSILICATE GLASS PROPERTIES |
| 5 | CRYSTAL CLEAR BENEFITS | 18 | CLOSED TUBULAR PHOTOBIO- REACTORS VERSUS OPEN PONDS |
| 7 | PRODUCT RANGE | 19 | BOROSILICATE GLASS VERSUS POLYMER MATERIALS |
| 8 | BOROSILICATE GLASS TUBING | 20 | PRESSURE DROP IN TUBULAR PHOTOBIOREACTORS |
| 10 | BOROSILICATE GLASS U- AND J-BENDS | 22 | REFERENCES OF GLASS TUBULAR PHOTOBIOREACTORS |
| 12 | COUPLINGS | | |
| 14 | BOROSILICATE GLASS MANIFOLDS | | |
| 15 | PACKAGING | | |

Algae Production Systems

Common photosynthetic algae cultivation systems are either open ponds or closed photobioreactors (PBRs).




Open ponds

Open ponds are typically built in circular or raceway configurations. The water is kept in motion, for example by paddle wheels. Open ponds are seemingly inexpensive and easy to build. However, poor light utilization, danger of contamination and high water evaporation are the main challenges, which lead to low biomass output per area and large water uptake. Some difficulties can be overcome by rooftops however this increases the costs further.

Closed system

Closed systems are dominated by tubular and flat-plate reactors. Other options are bags, coils or domes. Flat plate systems have received a lot of attention due to their large illuminated surface area, but the technology suffers from heating problems and a strong tendency to build up biofilm formations on the inner walls. Tubular systems on the other hand reduce these drawbacks while maintaining the advantages of optimal light input and high productivity.

Therefore closed tubular glass Photobioreactors (PBRs) with long lifetimes and easy cleanability, are very well suited for the highly reproducible cultivation of algae resulting in the highest possible growth rates. As such, tubular glass PBRs are best suited to provide biosecurity for high quality inoculum used in open ponds.



Upper photo
Haematococcus pluvialis, green phase
Lower photo
Haematococcus pluvialis, red phase

Right
Closed vertical tubular glass PBR system
© Jongerius ecoduna GmbH

Crystal Clear Benefits of Closed Tubular Glass PBR Systems



Biosecure ✓

protection against bio-contamination and culture crashes



Productive ✓

80 – 160 l/m² photoactive volume, PBR height up to 6 m



Cost efficient ✓

little maintenance and low total cost of ownership



Durable ✓

sustainable light transmission
T > 95 % (air – glass – water),
lifetime of 50 years and more



Resistant ✓

against chemicals, corrosion, sagging, scratches, UV-light




Food safe ✓

food and pharma grade




Product Range Overview

Helical System



- 1 U- or J-Bend
- 2 Coupling
- 3 Tubing

Fence System










horizontal
or vertical
orientation

- 1 Manifold
- 2 Coupling
- 3 Tubing


Product Range

Borosilicate Glass Tubing DURAN®

| Item No. | Joint Outside Diameter | Joint Wall Thickness | Tube Length | Volume per tube | Package Type* | Package Content | |
|----------|--|--|---|--|---|--|--|
| |  mm in |  mm in |  m ft |  l gal |  |  Number of Tubes |  Weight approx. kg approx. lb |
| 1535285 | 54 ± 0.54 2.13 ± 0.02 | 1.8 ± 0.12 0.07 ± 0.01 | 1.4 4.6 | 2.79 0.74 | carton | 9 | 8.3 18.3 |
| | | | | | pallet | 180 | 166 366 |
| 1522883 | | | 2.5 8.2 | 4.99 1.32 | carton | 12 | 19.7 43.5 |
| | | | | | pallet | 144 | 236.4 522 |
| 1523124 | | | 5.5 18 | 10.97 2.90 | wooden box | 56 | 202.6 430.8 |
| 1534297 | | | | | pallet | 238 | 861.2 1898.6 |
| 1500383 | 65 ± 0.65 2.56 ± 0.03 | 2.2 ± 0.18 0.09 ± 0.01 | 1.4 4.6 | 4.04 1.07 | carton | 9 | 12.2 26.9 |
| | | | | | pallet | 180 | 244 538 |
| 1511901 | | | 2.5 8.2 | 7.21 1.90 | carton | 9 | 21.8 48.0 |
| | | | | | pallet | 108 | 261.6 576 |
| 1459938 | | | 5.5 18 | 15.86 4.19 | wooden box | 36 | 191.6 422.3 |
| 1534302 | | | | | pallet | 165 | 877.9 1935.5 |

Other dimensions upon request.


*for explanation regarding package type please see *page 15*



Product Range

Borosilicate Glass U- and J-Bends DURAN®

Helical System



U-Bend



U-Bend long



J-Bend









Bend types are used for the following applications

- **U-Bend:** for gaps between tubes of typically 65 mm/2.56 in or more
- **U-Bend long:** used in combination with J- and U-Bends
- **J-Bend:** for smaller gaps down to 40 mm/1.57 in (couplings are shifted in every second row)










U-Bend

| Item No. | Joint Outside Diameter | Joint Wall Thickness | Joint U-Bend Width | Approx. U-Bend Height | Straight Side Length | Volume per Bend (approx.) | Package Type | Package Content |
|----------|---|---|---|---|--|---|---|--|
| |  mm in |  mm in |  mm in |  mm in |  mm in |  l gal |  |  Number of Tubes |
| 1534644 | 54 ± 0.65 2.13 ± 0.03 | 2.5 ± 0.14 0.10 ± 0.01 | 234 ± 2.0 9.21 ± 0.08 | 200 7.87 | > 45 > 1.77 | 0.85 0.22 | carton | 33 |
| | | | | | | | pallet | 396 |
| 1436672 | 65 ± 1.00 2.56 ± 0.04 | 2.8 ± 0.3 0.11 ± 0.01 | 245 ± 2.0 9.65 ± 0.08 | 200 7.87 | > 45 > 1.77 | 1.2 0.32 | carton | 21 |
| | | | | | | | pallet | 252 |

U-Bend long

| Item No. | Joint Outside Diameter | Joint Wall Thickness | Joint U-Bend Width | Approx. U-Bend Height | Straight Side Length | Volume per Bend (approx.) | Package Type | Package Content |
|----------|---|---|---|---|--|---|---|--|
| |  mm in |  mm in |  mm in |  mm in |  mm in |  l gal |  |  Number of Tubes |
| 1551070 | 65 ± 0.65 2.56 ± 0.03 | 2.8 ± 0.2 0.11 ± 0.01 | 245 ± 2.0 9.65 ± 0.08 | 300 11.81 | > 145 > 5.71 | 1.8 0.48 | carton | 21 |
| | | | | | | | pallet | 189 |

J-Bend

| Item No. | Joint Outside Diameter | Joint Wall Thickness | Joint U-Bend Width | Approx. U-Bend Height | Straight Side Length short leg | long leg | Volume per Bend (approx.) | Package Type | Package Content |
|----------|---|---|---|---|---|---|---|---|--|
| |  mm in |  mm in |  mm in |  mm in |  mm in |  mm in |  l gal |  |  Number of Tubes |
| 1657683 | 65 ± 0.65 2.56 ± 0.03 | 2.8 ± 0.2 0.11 ± 0.01 | 245 ± 2.0 9.65 ± 0.08 | 200 7.87 | > 45 > 1.77 | > 145 > 5.71 | 1.5 0.4 | carton | 21 |
| | | | | | | | | pallet | 189 |

Product Range

Couplings



SCHOTT tool kit



Open Standard Coupling




Open Standard Slim Coupling



Specially developed for tubular photobioreactors:
The couplings are designed for SCHOTT glass tubes with plain tube ends according to the product range shown in this brochure.

- Successfully tested for 10 years lifetime* regarding
 - 3 bar pressure resistance
 - UV-resistance
 - Regular cleaning cycles with various chemicals
- Fast installation allowing for reduced built up time of the reactor
- Easy to disassemble and re-use allowing for fast modification or extension of a reactor system
- Easy handling with pre-assembled devices and a special tool kit including a torque wrench (**Standard couplings: 40 Nm, Slim couplings 25 Nm**)
- Food grade


*see page 23




1 Standard, length 80 mm
Partition wall to assure glass separation and smooth transition (torque wrench 40 Nm needed)



2 Standard Slim, length 45 mm
(torque wrench 25 Nm needed)



3 Maintenance, length 80 mm
Allows easy exchange of tubes, no partition wall



4 Adapter, length 80 mm
Allows connection to periphery tubes with 2.5 in outer diameter

| Item No. | Outside Diameter Tube on Side 1 | Outside Diameter Tube on Side 2 | Package | Description | Package |
|----------|---|---|--|------------------------------|--|
| |  mm in |  mm in |  Number of Couplings | |  Weight approx. kg approx. lb |
| 1530116 | 54 2.13 | 54 2.13 | 24 bag | Standard ⁽¹⁾ | 6.9 15.3 |
| 1581056 | | | 960 pallet | | 277.7 612.2 |
| | 54 2.13 | 54 2.13 | 24 bag | Standard Slim ⁽²⁾ | 5.6 15.3 |
| | | | 1824 pallet | | 224.6 612.2 |
| 1463260 | 65 2.56 | 65 2.56 | 24 bag | Standard ⁽¹⁾ | 6.1 13.4 |
| | | | 960 pallet | | 219.4 483.7 |
| 1581035 | 65 2.56 | 65 2.56 | 24 bag | Standard Slim ⁽²⁾ | 5.9 13.4 |
| | | | 1824 pallet | | 236.2 483.7 |





Toolbox 54
Standard with torque wrench 40 Nm for closing of coupling and tools for opening

Toolbox 54
Slim with torque wrench 25 Nm for closing of coupling and tools for opening

Toolbox 65
Standard with torque wrench 40 Nm for closing of coupling and tools for opening

Toolbox 65
Slim with torque wrench 25 Nm for closing of coupling and tools for opening

Additional equipment

| Item No. | Outside Diameter Tube on Side 1 | Outside Diameter Tube on Side 2 | Package | Description | Package |
|----------|--|--|---|---|--------------------------------------|
| |  mm in |  mm in |  Number of Couplings |  | Weight approx. kg approx. lb |
| 1530120 | 54 2.13 | 54 2.13 | 4 carton | Maintenance kit (no partition wall) ⁽³⁾ | 0.9 2.0 |
| | | | 576 pallet | | |
| 1530105 | 65 2.56 | 65 2.56 | 4 carton | Maintenance kit (no partition wall) ⁽³⁾ | 1.0 2.2 |
| | | | 576 pallet | | |
| 1534828 | 65 2.56 | 63.5 2.50 | 4 carton | Adapter ⁽⁴⁾ | 1.0 2.2 |
| | | | 576 pallet | | |

Product Range

Borosilicate Glass Manifolds DURAN®

Fence System



Manifolds and couplings, fence system © Jongerius ecoduna GmbH



Manifold equipped with a flange



Manifold equipped with closed end

Manifolds are placed at the tops or at the ends of tubular PBR fences and function as U-Bends and in- and outlets.

- Biosecure and food safe, full glass solution
- Available with closed ends or with flange
- Outside diameter of arms: 54 mm or 65 mm for use with standard couplings
- Number of arms, distance between arms, total length etc. are customized with a minimum order quantity of 25 pieces

| Description* | Package | Package Content | Minimum Order Quantity (MOQ) |
|--------------------------------|---------|-----------------|------------------------------|
| Manifolds 10 arms closed 54 mm | | | 25 |
| Manifolds 5 arms flange 54 mm | | | 25 |
| Manifolds 8 arms closed 65 mm | | | 25 |
| Manifolds 4 arms flange 65 mm | | | 25 |
| Manifolds 6 arms open 65 mm | carton | 2 | 1 pallet |
| | pallet | 36 | |

*closed: both ends closed/flange: one side closed, one side flange/open: both sides flange

Packaging



Bags

- Couplings



Cartons

- Tubes, up to 2.5 m length
- U-Bends
- Couplings
- Manifolds



Wooden boxes

- Tubes, 5.5 m length
- smaller quantities



Special pallets

- Tubes, 5.5 m length

Borosilicate Glass

Properties

DURAN®


| | Metric | US |
|--|---|--|
| Coefficient of mean linear thermal expansion α acc. to DIN ISO 7991 | $3.3 \cdot 10^{-6} \text{ K}^{-1}$ (20 °C; 300 °C) | $3.3 \cdot 10^{-6} \text{ K}^{-1}$ (68 °F; 572 °F) |
| Transformation temperature T_g | 525 °C | 977 °F |
| Density ρ at 25 °C | 2.23 g · cm ⁻³ | 139.2 lb · ft ⁻³ |
| Modulus of elasticity E (Young's modulus) | $63 \cdot 10^3 \text{ N} \cdot \text{mm}^{-2}$ | $91 \cdot 10^5 \text{ lb} \cdot \text{in}^{-2}$ (psi) |
| Poisson's ratio μ | 0.20 | 0.20 |
| Thermal conductivity λ_w at 90 °C | $1.2 \text{ W} \cdot \text{m}^{-1} \cdot \text{K}^{-1}$ | $0.69 \text{ Btu} \cdot \text{hr}^{-1} \cdot \text{ft}^{-1} \cdot ^\circ\text{F}^{-1}$ |
| Refractive index ($\lambda = 587.6 \text{ nm}$) n_d | 1.473 | 1.473 |
| Stress-optical coefficient (DIN 52 314) K | $4.0 \cdot 10^{-6} \text{ mm}^2 \cdot \text{N}^{-1}$ | $4.0 \cdot 10^{-6} \text{ mm}^2 \cdot \text{N}^{-1}$ |

| Chemical composition | | | |
|----------------------|-------------------------------|---|--------------------------------|
| SiO ₂ | B ₂ O ₃ | Na ₂ O + K ₂ O | Al ₂ O ₃ |
| 81 | 13 | 4 | 2 |

main components in approx. weight %

| Chemical resistance | |
|--------------------------------|-----------|
| Hydrolytic Class (DIN ISO 719) | HGB 1 |
| Acid Class (DIN 12116) | Class S 1 |
| Alkali Class (DIN ISO 695) | Class A 2 |

Transmission



Transmission of DURAN® glass (d=2.2mm) in configuration air/glass/air.

Note

When the glass tube is filled with water, the transmission increases from about 92 % to 95.6 % due to reduced reflection losses at the inner glass/water interface.

Pressure resistance of tubing made of borosilicate glass

The following formulas apply to stress free, pristine tubing and cylindrical hollow bodies with a circular profile, uniform wall thickness with open ends, free from thermal load, under internal positive pressure.

Estimation of the maximum pressure resistance (p)

$$p = \frac{WT \cdot 140 \text{ bar}}{OD - WT}$$

OD = Outside Diameter in [mm]
WT = Wall Thickness in [mm]
p = Pressure Resistance in [bar]

Estimation of the minimum wall thickness (WT)

$$WT = \frac{OD \cdot p}{140 \text{ bar} + p}$$

$$\frac{K}{S} = 70 \text{ bar}$$

Permissible load referring to standard DIN EN 1595: "Pressure equipment made from borosilicate glass 3.3 – general rules for design, manufacture and testing"

Other points to be considered:

- AD 2000-leaflet N 4, edition 2000-10: Pressure vessels made of glass, with encl. 1, edition 2000-10: evaluation of faults in walls of glass pressure containers
- AD 2000-leaflet B 1, edition 2000-10: cylindrical and spherical shells under internal pressure overload

According to DIN EN 1595: "Pressure equipment made from borosilicate glass 3.3 – general rules for design, manufacture and testing", DURAN® approved material and may be used for the construction of pressure equipment.

Features and Benefits

Closed Tubular Photobioreactors versus Open Ponds



Closed tubular photobioreactor



Open pond

Contamination

- Very low risk of contamination compared to open ponds, where other microorganisms or insects have easy access
- No limitation regarding the algae species that can be grown, also due to effective blocking against competing organisms

Productivity

- Higher productivity in terms of mass per area and day
- Significantly higher volumetric productivity

Algae concentration at harvest

- Notably higher concentration in terms of mass per liter
- More efficient harvesting procedure

Water loss

- No evaporation within closed system compared to open ponds, which can lose significant water amounts, resulting in salinization
- Water loss is limited to external factors, such as the cooling process

Biomass quality

- Biomass quality is highly reproducible due to excellent process control of tubular PBR systems
- High value products or high quality inoculum can be produced with optimum reliability

Production flexibility

- Easy cleanability allowing for defined initial status any time, thus switching algae species is possible and secure


Use of GMO* for improved production process

- GMO production is possible with closed reactor design


*GMO = Genetically Modified Organism

Borosilicate Glass versus Polymer Materials


| | |
|---------------------------|---|
| Light transmission | <ul style="list-style-type: none"> ■ Excellent light transmission (see page 17 for details) ■ No solarization or browning effect ■ No UV-protective additive or coating necessary to secure material properties |
| Fire protection | <ul style="list-style-type: none"> ■ Glass does not burn or give off toxic fumes |
| Leaching | <ul style="list-style-type: none"> ■ Glass is a chemically highly resistant material. With plastic tubing, depending on the polymer type, monomers or oligomers of hazardous substances such as Bisphenol molecules can be leached into the algae culture. |
| Cleaning | <ul style="list-style-type: none"> ■ Mechanical stability allows continuous in-line cleaning with polymer pellets ■ Chemical stability allows cleaning in place (CIP) ■ Lower material and maintenance costs compared to quality polymer tubes |
| Thermal stability | <ul style="list-style-type: none"> ■ No need for expansion loops due to low thermal expansion Example: For 5.5 m long tubes and a temperature increase of 20 °C/36 °F the expansion of Borosilicate glass is only 0.36 mm/0.01" while polymers expand from 3.3-8.8 mm/0.13"-0.35" depending on polymer type. |
| Cost saving | <ul style="list-style-type: none"> ■ Glass components can last fifty years and longer ■ Reduced number of rack poles. High mechanical stability allows increased support distances without significant sagging of tubes (see picture on right) |
| Sagging | <ul style="list-style-type: none"> ■ No permanent deformation of glass tubes in contrast to polymer tubes |



Borosilicate glass



Polymer




Sagging of water filled tubes

(outer diameter 65 mm, wall thickness 2.2 mm, length 2.75 m).


The sag of the glass and polymer tubes is 0.5 mm and 6.6 mm, respectively. The polymer tube would need to be supported every 1.5 m for the same sag as the glass tube.

Pressure Drop in Tubular Photobioreactors

For the optimum construction of a photobioreactor the expected pressure loss in the system must be known. This information is relevant for the lay-out of the ideal pump design. For the calculation of the pressure drop of the entire photobioreactor system the pressure drops of the individual components can be added together. The individual pressure drops must be calculated for the desired velocity, u . In the following table (p. 21) the individual pressure drops of representative glass components are shown for a typical flow velocity of 0.7 m/s. The pressure drops depend on the dimensionless zeta-values, which slightly decrease at larger velocities. Please contact SCHOTT technical service for assistance and additional simulations.



Dean-vortex appearance in an U-Bend – computer simulation (ANSYS® CFX® 14.5.7)



Pressure loss

In general, the pressure drop can be calculated for any velocity using the following formula.

$$\Delta p = \zeta \cdot \frac{\rho}{2} \cdot u^2$$

Δp : pressure loss
 ζ : pressure loss number (zeta)
 ρ : algae culture density
 u : linear velocity of algae culture

| $u = 0.7 \text{ m/s}$ | ζ | $\Delta p \text{ [Pa]}$ |
|---|---------|-------------------------|
| Round Tube ($D = 65 \text{ mm}$, $WT = 2.2 \text{ mm}$, $L = 5.5 \text{ m}$) | 1.96 | 480 |
| U-Bend ($D = 65 \text{ mm}$, $WT = 2.8 \text{ mm}$) | 0.252 | 62 |

Pressure drops of a tube and a U-Bend at the given velocity, u . D is the outer diameter, WT the wall thickness. The algae culture's density was approximated with $\rho = 1 \text{ g/cm}^3$.

Electrical power

The electrical power of the pumps, P_{el} , scales with the pressure drop and the volume flow, Q :

$$P_{el} = \frac{\Delta p \cdot Q}{\eta_p}$$

P_{el} : electrical power
 Δp : sum of pressure loss in Pa
 Q : volume flow rate in m^3/s
 η_p : pump efficiency at operating point ($0 < \eta_p < 1$)

Research done in cooperation with: Institute of Fluid Dynamics, LSTM, Technical Faculty, University Erlangen-Nuremberg, Germany

References

of Glass Tubular Photobioreactors

SCHOTT has formed alliances and partnerships all over the world. This allows us to provide complete tubular photobioreactors according to your needs. Please contact us for further details.

1 © Algalíf Iceland ehf., Iceland


2 © Varicon Aqua Solutions Ltd., UK, Phyco-Flow PBR at OP Bio, Japan


3 © A4F-Algae for Future, Portugal

4 © Algatechnologies Ltd., Israel

5 © Clearas Water Recovery Inc., USA

6 © Jongerius ecoduna GmbH, Austria









Technical Terms of Supply

Detailed information on permissible faults, definition of faults, testing methods and testing units are available upon request. Reduced tolerances are also available upon request. Regarding quality issues the relevant “Technical Terms of Supply” for the application apply to all sales and are binding unless separate written agreements with respect to specification have been agreed upon.

SCHOTT® is a registered trademark of SCHOTT.
DURAN® is a registered trademark of DWK Life Sciences GmbH.

We thank our customers and partners for their kind assistance in providing product samples and photos.

SCHOTT declines any liability with respect to the correctness or completeness of any information or data that is expressly or implicitly contained in this brochure.

| | |
|---|---|
|  <p>Management System ISO 9001:2015 ISO 15378:2017 www.tuv.com ID 9108654402</p> <p>Germany: Site Mitterteich + Mainz</p> |  <p>Management System ISO 9001:2015 ISO 15378:2017 www.tuv.com ID 9105082675</p> <p>India: Site Baroda</p> |
|  <p>Management System ISO 9001:2015 ISO 15378:2017 www.tuv.com ID 9105082675</p> <p>Brazil: Site Rio de Janeiro</p> |  <p>Management System ISO 9001:2015 ISO 15378:2017 www.tuv.com ID 9000019772</p> <p>China: Site Jinyun</p> |

schott.com/pbr

SCHOTT North America, Inc., 2 International Drive, Suite 105
Rye Brook, NY 10573, USA