

SCHOTT TopLyo®

Datasheet



SCHOTT TopLyo® – Pharmaceutical glass vial with hydrophobic coating to avoid fogging

- For stability reasons, more than 50 % of all biologic drugs are lyophilized.
- Fogging is a widely known undesirable phenomenon that occurs during lyophilization, which results in elevated levels of rejects.
- Overfilling vials to compensate for drug loss through unsatisfactory residual emptying after reconstitution leads to higher costs.
- SCHOTT TopLyo® is unique in combining hydrophobic behavior and the avoidance of free silicone.
- Inner coating is applied using patented and proven plasma impulse chemical vapor deposition (PICVD) technology.
- An additional “release criterion” has been specifically developed.



Si-O-C-H layer
applied via PICVD



Layer thickness of ~ 40 nm



Long-term
layer stability



Coating bonds covalently
to the glass substrate



Suitable for
depyrogenation



Dense coating
(i.e. non-porous)



Contact angle for water
> 90° (hydrophobic surface
without silicone)



Stable after washing
process



Compliant with all current
standards, such as Ph.
Eur, USP, JP, and CP

No fogging

Particularly suitable for antibody drug conjugates (ADCs) thanks to hydrophobic inner surface

Elegant lyo cake

Reduced rejects due to improved lyo cake aesthetics

Improved emptying


Less residual volume so no overfilling necessary

Drug stability


No free silicone thanks to residual free technology

SCHOTT
PHARMA

Release test: Drain-off test for hydrophobicity




Verifications: Improved emptying



Result
Less residual volume after reconstitution with SCHOTT TopLyco vials

Verifications: No fogging and elegant lyo cake



Method

10R vials Type I glass and SCHOTT TopLyco vials.
5.0 ml formulation dried in 10R vials with different surfaces.
0.15 mg/ml human growth hormone, 40 mg/ml mannitol, and 10 mg/ml sucrose. Phosphate/glycine buffer (pH 7.0). Sterilization using 0.2 µm PES Filter, 25 °C, 30 min.

Result

Less cake disruption and dry material pulling from the edge with SCHOTT TopLyco vials.

Verifications: Stress tests have proven stability

	TopLyco® 10R	TopLyco® 15R	TopLyco® 10R, depyrogenated	TopLyco® 15R, depyrogenated
ca. 5 mm*	average 103	102	99	101
	stand. dev./ range ± 2 / ± 4	± 2 / ± 4	± 2 / ± 5	± 2 / ± 4
ca. 15 mm*	average 102	102	98	100
	stand. dev./ range ± 2 / ± 5	± 2 / ± 4	± 2 / ± 5	± 1 / ± 3
ca. 25 mm*	average 106	103	101	101
	stand. dev./ range ± 4 / ± 9	± 1 / ± 3	± 3 / ± 5	± 1 / ± 3

Method

TopLyco vials: 10R (> four years of storage) and 15R (three months storage).
Contact angle measurement at three lateral positions (bottom, middle, and neck area)*:
Reference vs. depyrogenated (30 min at 330 °C). 15 vials measured per sample type.

Result

All analyzed vials show hydrophobic behavior with stable contact angle > 90°. No significant differences were observed for different storage times.

General ordering information

Quality level	TopLine with additional release test										
Packaging	<ul style="list-style-type: none"> Tray with optional divider Pre-washed and pre-sterilized: adaptiQ® (tray, cup nest, clip nest*) 										
Palletizing	Standard Euro pallet (1200 x 800 mm) contains 15-27 layers of nine trays each										
Formats	2R	3R	4R	6R	8R	10R	15R	20R	25R	30R	50R
Pieces per tray	344	344	344	186	186	154	154	95	95	95	40

*50R not available in clip nest configuration

Many configurations are available in small quantities as "Fast Track Articles".

Visit our online shop or speak to your sales representative for more information.

FIND OUT MORE:
www.schott-pharma.com/vials

**SCHOTT
PHARMA**