

Non-Invasive In-Situ Fluorescence Spectroscopy: Flexible Monitoring of Complex Bioprocesses in Stirred Tank Fermenters


Christian Ott¹, Dr. Ivan Schlembach², Kristina Gruber¹, Dr. Lars Regestein²

¹ SCHOTT AG, Landshut, Germany | ² Leibniz Institute for Natural Product Research and Infection Biology – Hans Knöll Institute, Jena, Germany

1. INTRODUCTION

- **Challenge:** With its efficiency-boosting potential, biocatalysis is becoming more important to pharma and biotech industries. Economically making the needed enzymes is difficult because insoluble media components hamper continuous process control.
- **Approach:** In-situ fluorescence spectroscopy can analyze biomass growth, substrate consumption and product formation during bioprocesses and is less effected by insoluble particles compared to optical density or scattered light measurements, which cannot distinguish between light absorption by cells and cellulose fiber. However, the practicality of fluorescence spectroscopy depends on the ability to mount and unmount the spectrometer probe during cultivation without compromising sterile integrity.
- **Enabler:** To overcome this challenge, SCHOTT developed the ViewPort®, which provides a sterile-safe interface that enables in-situ bioprocess monitoring through a hermetically sealed optical window.

2. PRINCIPAL SETUP and ASSEMBLY



3. USE CASE: IN-SITU MONITORING of ENZYME PRODUCTION in a SUSPENSION

- We used a stirred-tank fermenter (2.5 L fermentation volume), ViewPort® optical interface (SCHOTT, Germany) and Carry Eclipse Fluorescence spectrometer (Agilent, USA) with fiber optic probe (Fig. 1).
- The fungus *Trichoderma reesei* was cultivated in a medium containing 5 g/L glucose, 30 g/L alpha cellulose as carbon sources and 2 g/L peptone.
- The fungus was tagged to express mCherry fluorescent protein as marker to follow growth.
- We compared in-situ fluorescence data with data from online off-gas monitoring (non-optical approach, real-time estimates of cell growth and substrate consumption).
- An excitation vs. emission scan (Fig. 2) was used to detect biomass autofluorescence signals.



Fig. 1: Experimental setup





Fig. 2: Spectrally resolved in-situ fluorescence measurement for discrimination

4. RESULTS



- Standard off-gas analysis parameter CTR correlated with CDM until 100 h of cultivation (exponential cell growth); in-situ mCherry fluorescence signal correlated until 144 h (Fig. 3A); difference is likely due to metabolism change.
- Recording autofluorescence signals for tryptophan and flavin; flavin correlated better with cell growth than mCherry or tryptophan (Fig. 3B); setup's sensitivity allowed detection of autofluorescence signals with low response (flavin).

5. SUMMARY



Autofluorescence of flavin is a suitable process parameter (vs. mCherry fluorescence)



The in-situ set-up is sensitive enough to detect autofluorescence signals



Autofluorescence detection enables marker-free bioprocess monitoring – a prerequisite for production scales



ViewPort® enables sterile-safe, sensitive and flexible in-situ monitoring

Fig. 3: Fluorescence-based characterization of growth and protein formation of *T. reesei* RUT-C30 mCherry

Link to TechnoPharm article (in German):

