




# EVERIC<sup>®</sup> Iyo

## Pharmaceutical glass vial with hydrophobic coating to avoid fogging



- For stability reasons, more than 50 % of all biologic drugs are lyophilized.
- Fogging is a widely known undesirable phenomenon that occurs during lyophilization, which results in elevated levels of rejects.
- Overfilling vials to compensate for drug loss through unsatisfactory residual emptying after reconstitution leads to higher costs.
- EVERIC<sup>®</sup> Iyo is unique in combining hydrophobic behavior and the avoidance of free silicone.
- Inner coating is applied using patented and proven plasma impulse chemical vapor deposition (PICVD) technology.
- An additional "release criterion" has been specifically developed.



Si-O-C-H layer  
applied via PICVD



Layer thickness  
of ~ 40 nm



Long-term  
layer stability



Coating bonds covalently  
to the glass substrate



Suitable for  
depyrogenation



Dense coating  
(i.e. non-porous)



Contact angle for water  
> 90° (hydrophobic  
surface without silicone)



Stable after  
washing process



Compliant with all current  
standards, such as Ph. Eur,  
USP, JP and CP



### No fogging

Particularly suitable for antibody-  
drug conjugates (ADCs) thanks to  
hydrophobic inner surface



### Elegant Iyo cake

Reduced rejects due to  
improved Iyo cake aesthetics



### Improved emptying

Less residual volume so  
no overfilling necessary




### Drug stability


No free silicone thanks to  
residual free technology

**SCHOTT**  
**PHARMA**

Release test: Drain-off test for hydrophobicity



Verifications: Improved emptying



Verifications: No fogging and elegant lyo cake



Verifications: Stress tests have proven stability

		EVERIC® lyo 10R	EVERIC® lyo 15R	EVERIC® lyo 10R, depyrogenated	EVERIC® lyo 15R, depyrogenated
ca. 5 mm*	average	103	102	99	101
	stand. dev./range	± 2 / ± 4	± 2 / ± 4	± 2 / ± 5	± 2 / ± 4
ca. 15 mm*	average	102	102	98	100
	stand. dev./range	± 2 / ± 5	± 2 / ± 4	± 2 / ± 5	± 1 / ± 3
ca. 25 mm*	average	106	103	101	101
	stand. dev./range	± 4 / ± 9	± 1 / ± 3	± 3 / ± 5	± 1 / ± 3

Method

10R vials Type I glass and EVERIC® lyo vials.  
5.0 ml formulation dried in 10R vials with different surfaces.  
0.15 mg/ml human growth hormone, 40 mg/ml mannitol, and 10 mg/ml sucrose. Phosphate/glycine buffer (pH 7.0).  
Sterilization using 0.2 µm PES Filter, 25°C, 30 min.

Result

Less cake disruption and dry material pulling from the edge with EVERIC® lyo vials.

Method

EVERIC® lyo vials: 10R vial (> four years of storage) and 15R vial (three months storage). Contact angle measurement at three lateral positions (bottom, middle, and neck area)\*: Reference vs. depyrogenated (30 min at 330 °C). 15 vials measured per sample type.

Result

All analyzed vials show hydrophobic behavior with stable contact angle > 90°. No significant differences were observed for different storage times.

General ordering information											
Quality level	TopLine with additional release test										
Packaging	<ul style="list-style-type: none"><li>Tray with optional divider</li><li>Pre-washed and pre-sterilized: adaptiQ® (tray, cup nest)</li></ul>										
Palletizing	Standard Euro pallet (1200 x 800 mm) contains 15–27 layers of nine trays each										
Formats	2R	3R	4R	6R	8R	10R	15R	20R	25R	30R	50R
Pieces per tray	344	344	344	186	186	154	154	95	95	95	40


Many configurations are available in small quantities as “Fast Track Articles”.

Visit our Online Shop or speak to your sales representative for more information.



# EVERIC<sup>®</sup> amber

Reliable light protection and regulatory compliance



- While every protein is sensitive to light, the more complex the molecule, the higher the possibility of light-induced aggregation or degradation.
- Antibody-drug conjugates (ADCs) are particularly sensitive to light due to their fragility. Even light exposure during filling prior to lyophilization can cause issues.
- Adapting the lighting set-up during filling is cumbersome and may not be sufficient.
- Existing solutions do not comply with USP, Ph. Eur., and JP light transmission requirements, so multiple primary packaging articles may be necessary.
- EVERIC<sup>®</sup> amber offers a solution that provides reliable light protection while remaining compliant with USP, Ph. Eur., and JP standards.



FIOLAX<sup>®</sup> amber glass tubing –  
absorption of light



Improved forming process: tight  
control of dimensional parameters  
and temperature ranges



Unchanged glass composition –  
Type I Borosilicate Glass



Full compliance with light  
transmission requirements of  
Ph. Eur., USP and JP



## Light protection

Protection against ultraviolet rays  
and short-wave visible light



## Global compliance

Suitable for pharmaceutical companies  
operating in multiple regulatory environments

**SCHOTT**  
PHARMA

Release test


01 Mechanical cut & cleaning

- Vial is:
- 1) Cut in vertical direction
  - 2) Cleaned to get a quasi-plan, contaminant-free surface




02 Spectral screening

Light transmission is measured through the middle of the vial body



03 Spectral evaluation

Certified release criteria acc. to wall thickness and filling volume



Verifications: The right production set-up is key to achieving global compliance

Annealing after hot-forming is necessary to avoid residual stress, but it influences transmission.

Light transmission compliance with Ph. Eur., USP and JP standards can be achieved for amber vials with a 1 mm thick wall (2-10R) via tight process control.




Figure 1: Light transmission compliance for 10R vial (wall thickness 1 mm)

By reducing the wall thickness of a 20R vial to 1 mm, compliance with USP, Ph. Eur., and JP standards can be achieved.

Strength-optimized parameters with tight dimensional tolerances, dedicated inspection and packaging with dividers ensure strength similar to conventional vials.




Figure 2: Strength tests for 20R vial with a reduced wall thickness (1 mm)

General ordering information					Order now: www.schott-pharma.com/vials
Quality level	StandardLine & release test				
Packaging	<ul style="list-style-type: none"><li>▪ Tray (with divider for 10R and 20 ml format)</li><li>▪ Pre-washed &amp; pre-sterilized: adaptiQ® (tray, cup, nest)</li></ul>				
Possible combinations	EVERIC® amber can be combined with EVERIC® Iyo				
Palletizing	A standard Euro Pallet (1200 x 800 mm) contains 15 – 27 layers of 9 trays each				
Formats	2R	6R	10R	20 ml	
Pieces per tray	344	186	104	77	