

FIOLAX® clear with Controlled Hydrolytic Resistance (CHR)

General Product Information

The increasing demand of a reduced sodium release from glass led to the development of a FIOLEX® clear tubing with controlled hydrolytic resistance (CHR) of the inner surface. This quality feature enables the converting process based on a known and defined low sodium release from the starting material. For the first time, not only the hydrolytic resistance of the glass grains, but also from the inner surface is monitored and specified for the glass tubing – therefore: The best quality right from the start! The composition, chemical as well as physical properties of Type I borosilicate FIOLEX® clear glass remain fully untouched, hence FIOLEX® clear CHR can be used without any need to change registration files.

Value-adding Product Benefits



Guaranteed quality right from the start

The specification of a limit for the hydrolytic resistance of the inner surface of the tubing sheds light into this so far unknown aspect – always guaranteed by a certificate.



Applicable to registered products

FIOLEX® clear tubing can be replaced by FIOLEX® clear CHR tubing without the need of a new registration of the pharmaceutical product.



Superior product quality especially for sensitive formulations

FIOLEX® clear CHR is especially beneficial for highly demanding formulations, which are sensitive towards extractables such as sodium, e.g. biotech drugs or pH sensitive / unbuffered formulations.



SCHOTT


FIOLAX® clear CHR

Facts at a Glance

Verification


The hydrolytic resistance of the inner surface of a container is composed by two different aspects – the initial level of the tubing as well as the influence of the converting process. Whereas many efforts were made optimizing the production process of a container, the level of the glass tubing at the beginning mostly remains an unknown aspect.

By implementing a tighter process control during tube production as well as the introduction of the inner surface test parallel to the tubing production, it is now possible to guarantee a tubing with a specified, low sodium release of the inner surface. Exemplarily, measurement data of tube sections acc. to ISO 4802-2:2016 of FIOLAX® clear CHR are shown in the graph (referring to 2R vial dimensions: outside diameter 16.00 mm, wall thickness 1.00 mm). The data reveals that it is possible to specify a limit of 0.70 µg/ml for this tubing dimension, which is below 25 % of the ISO 4802-2 limit for a 2R vial.



- Typical range for a 2R vial
- Typical range for 2R vial tubing
- Limit ISO 4802-2 for final 2R vial
- SCHOTT Limit FIOLAX® clear CHR 2R vial tubing
- Exemplary measurement results of FIOLAX® clear CHR 2R vial tubing

Availability





FIOLAX® clear CHR tubing is produced according to the current Technical Performance Specification (TPS) for pharmaceutical Type I glass tubing including SCHOTT perfeXion® technology. It is available within our **Best value**  customizing options for FIOLAX® clear CHR for an outside diameter range of 6.85 mm up to 50.00 mm.

Test and Specification

For FIOLEX® clear CHR a limit for the sodium extraction from the inner surface acc. to ISO 4802-2:2016 is specified. Since the release depends on the tubing diameter, the limit is specific for each dimension and is adapted to tubing. The limits of standard sizes are displayed in the table below, customized dimensions are also available upon request. After autoclaving with purified water for 1 h at 121 °C the amount of sodium oxide in the extraction solution is determined by means of flame spectrometry and expressed as micrograms of Na₂O per milliliter of extraction solution [µg/ml], which is equivalent to parts per million [ppm].



For each delivery of FIOLEX® clear CHR tubing a Certificate of Conformity (CoC) will be supplied.

	Container Size [ml]*	Outside Diameter [mm]	Inside Diameter [mm]	Wall Thickness [mm]	Limit Value for HC _F 1 acc. to ISO 4802-2 [µg/ml]	Limit Value for FIOLEX® clear CHR Tubing [µg/ml]
 Syringe	0.50	6.85	4.65	1.10	5.00	2.20
	1.00	8.15	6.35	0.90	4.50	1.60
	1.00 / 2.00 / 2.25 / 3.00	10.85	8.65	1.10	4.50 / 4.10 / 4.10 / 3.20	1.20
	5.00	14.45	11.85	1.30	2.50	0.90
	10.00	17.05	14.25	1.40	2.00	0.80
	20.00	22.05	19.05	1.50	1.50	0.60
 Cartridge	1.50	8.65	6.85	0.90	4.50	1.40
	3.00	11.60	9.70	0.95	3.20	1.00
	10.00	22.05	19.05	1.50	2.00	0.60
 Vial	2R / 3R / 4R	16.00	14.00	1.00	3.20 / 3.20 / 2.50	0.70
	6R / 8R	22.00	20.00	1.00	2.50 / 2.00	0.50
	10R / 15R	24.00	22.00	1.00	2.00 / 2.00	0.50
	20R / 25R / 30R	30.00	27.60	1.20	1.50 / 1.50 / 1.50	0.50
	50R	40.00	37.00	1.50	1.50	0.50
	100R	47.00	43.60	1.70	1.20	0.50
 Ampoule	1.00	10.75	9.75	0.50	4.50	1.00
	2.00	10.75	9.75	0.50	4.10	1.00
	5.00	14.75	13.65	0.55	2.50	0.80
	20.00	22.50	21.10	0.70	1.50	0.50

* acc. to ISO 11040-4 (syringes), 8362-1 (vials), 9187-1 (ampoules)

All specifications are subject to change without prior notice. This datasheet or any extracts thereof may only be used in other publications with express permission of SCHOTT.

[schott.com](https://www.schott.com)

SCHOTT AG, Erich-Schott-Strasse 14, 95666 Mitterteich, Germany

SCHOTT