
EasyLED User Manual

DE | EN | FR | ES | IT

SCHOTT
glass made of ideas

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

Device overview

- ❶ Electronic light intensity control
- ❷ Mains connection socket
- ❸ Indicator light
- ❹ On/off switch
- ❺ Clamping screws
- ❻ Jog dial (Only with EasyLED Ring Light Plus)
- ❼ Information label

- ❶ Electronic light intensity control
- ❷ On/off switch
- ❸ LED head
- ❹ Mains connection socket
- ❺ Information label
- ❻ Heat sink
- ❼ Indicator light
- ❽ Attachment adapter socket

Inhalt

1. Important information	16
2. Operation Ring Lights	18
2.1 Assembly	18
2.2 Power Connection	18
2.3 Start-up procedure	18
2.4 Segment setting	18
2.5 Single step and continuous rotation	19
3. Operation Spot Lights	20
3.1 Assembly	20
3.2 Power connection	20
3.3 Start-up procedure	21
3.4 Light intensity setting	21
3.5 Operating mode selection	21
4. Transmitted Light Stage	22
4.1 Assembly	22
4.2 Power connection	22
4.3 Start-up procedure	22
5. Maintenance	23
6. Troubleshooting	23
7. Technical data	24

1. Important information

Symbols used

Symbol	Meaning
	Warning of danger (caution, obey documentation)

Intended use

The EasyLED series illumination system is intended for use in industrial and laboratory applications¹. The system is used to illuminate objects in stereo-microscopy. The Transmitted Light Stage is designed to operate in a microscope stand. Optional adapter rings can be used, various sizes available. Only using SCHOTT accessories will guarantee proper function, safety and optimum light output.

Safety information

Please read and follow these instructions carefully. Ensure that every operator of the device has quick access to these operating instructions.

The instrument's safety cannot be guaranteed if the instructions are not followed. The manufacturer is not liable for any damage caused by non-observance of these operating instructions.

- Operate the lighting system only in dry rooms and not in areas subject to explosive hazards.
- Operate the lighting system only with the voltage specified on the information label.
- Never look directly into the light emission aperture of the lighting system during operation (risk of eye injury)! According to the EN 62471:2008 standard, the LED lights in the EasyLED series are classified as Risk Group 0 (transmitted light), Risk Group 1 – low risk (ring lights), and Risk Group 2 – medium risk (spot lights).
- Never cover the lighting system when system is turned on (fire hazard).
- Never cover the light exit opening of the illumination system with your hand or other parts of your body.

¹ You can obtain an individual spectrum of your light source at info.microscopy@schott.com

- When illuminating heat-sensitive or flammable light-absorbing objects, pay particular attention to selecting the distance of the light exit or the lamp brightness so that no thermal damage occurs to the object. Since light-absorbing materials have the physical property of converting incident light into heat, damage can occur to heat-sensitive or flammable light-absorbing materials.
- Make sure that all unused light heads are always at a safe distance – at least 10 cm – from heat-sensitive or flammable light-absorbing materials when the Spot Light is switched on (avoid possible fire hazard). Ensure the above safe distance e.g. for dark/colored textiles and dark/colored wood or plastic surfaces.
- Reduce the brightness and duration of illumination of biological tissue to that which is strictly necessary to avoid unnecessary exposure of biological tissue to visible light illumination.
- Safe disconnection from the power supply is only achieved by pulling out the plug.
- Do not open or disassemble the devices. Refrain from making technical modifications to the devices. Repairs may only be carried out by the manufacturer or customer service centers authorized by the manufacturer.

2. Operation Ring Lights

2.1 Assembly

Attach the Ring Light to the objective of the microscope using the two clamping screws so that the green indicator light faces the user. The inner diameter of the Ring Light is 66 mm. Appropriate adapter for other diameters are available for connection to objectives.

2.2 Power connection

Insert the hollow plug of the enclosed power supply unit into the connection socket of the Ring Light provided for this purpose. Make sure that your Ring Light is only operated in the specified voltage range and with the enclosed power supply unit.

2.3 Start-up procedure

You can switch the Ring Light on and off by pressing the main switch. If the device is switched on, the green control lamp lights up. Set the desired brightness on the electronic light intensity control. To disconnect the power supply, pull out the mains plug.


2.4 Segment setting

The EasyLED Ring Light Plus can be controlled segment by segment.

When the segment setting is matched to the object being viewed, object contrast can be achieved that far exceeds the contrast when illuminated with a non-segmented ring light.

The different segment modes are set by pushing the electronic jog dial. Repeated pushing of the jog dial changes the segment mode of the LED unit.

„Full circle -> semicircle -> quarter circle -> double quarter circle -> full circle“.
The Ring Light starts in full-circle mode after each switch-off.

With the Ring Light Plus, you can rotate the segment setting step by step or set it to continuous rotation. Rotation mode increases the impression of three-dimensionality of structured objects, especially when viewed on the monitor.

By rocking the jog dial clockwise in the segment setting, the set segments move correspondingly clockwise in 1/8 circle steps. Rocking the jog dial counterclockwise in the segment setting causes the set segments to move counterclockwise accordingly.

Prolonged dragging of the jog dial to the right or left activates the continuous rotation of the LED unit in the corresponding direction.

The automatic rotation stops as soon as the jog dial is switched arbitrarily.



2.5 Single step and continuous rotation

3. Operation Spot Lights

3.1 Assembly

The EasyLED Spot Light Plus can be mounted on all common microscope or swing-arm stands. Adapters are included with the light source or are available from the dealers or manufacturer.

Insert the mounting adapter (A) into the receptacle on the underside of the light source as far as it will go and tighten the two mounting screws (B) using an Allen key (3 mm). You can now attach the Spot Light to the microscope stand or swivel arm stand by clamping it using the mounting adapter. Keep in mind, you may have to remove the microscope body or the carrier beforehand and mount it again afterwards.

3.2 Power connection

Insert the hollow plug of the enclosed power supply unit into the socket provided for this purpose on the control unit. **Make sure that your Spot Light is only operated in the specified voltage range and with the enclosed power supply unit.**

You can switch the Spot Light on and off by pressing the main switch. If the device is switched on, the blue control lamp lights up. To disconnect the power supply, pull out the plug!.

Adjust the light intensity by turning the light intensity control. The scale is divided into six brightness levels, from 0 to 5. The highest light intensity is achieved at full deflection in position 5.

Attention: In the control position 0 the device is not switched off, therefore the blue control light is on. To switch off the Spot Light, press the main switch.

By repeatedly pressing the main switch, you can select between the following operating modes for the Spot Light Plus Double:

Light head 1 ON, light head 2 ON, both light heads ON, both light heads OFF, etc.

3.3 Start-up procedure

3.4 Light intensity setting

3.5 Operating mode selection

4. Transmitted Light Stage

4.1 Assembly

Place or attach the Transmitted Light Stage in the recess provided on the microscope stand so that the green indicator light faces the user. The basic version of the Transmitted Light Stage has a diameter of 84 mm. Adapters are available for other commonly used diameters to accommodate different microscope brands.

4.2 Power connection

Insert the hollow plug of the enclosed power supply unit into the socket provided for this purpose on the control unit. Make sure that your Transmitted Light Stage is only operated in the specified voltage range and with the enclosed power supply unit.

4.3 Start-up procedure

You can switch the Transmitted Light Stage on and off by pressing the main switch. If the unit is switched on, the green control lamp lights up. To disconnect the power supply, pull out the power plug.

The EasyLED series system is maintenance-free.

Use a soft dry cloth or commercially available plastic cleaning cloths for external cleaning of the device.

The use of cleaning agents, alcohol or any other chemicals is not recommended.

If you are unable to remedy the fault using the measures described below, please contact your specialist dealer or the nearest SCHOTT agency.

5. Maintenance

6. Troubleshooting


Typ of fault	Possible cause	Remedy
1. Lamp out	<ul style="list-style-type: none"> • Instrument not switched on • Plug not in socket • Power not provided to instrument 	<ul style="list-style-type: none"> • Switch instrument on • Plug in socket • Check instrument is connected to power
2. Light intensity dropping	<ul style="list-style-type: none"> • Lighting system overheated 	<ul style="list-style-type: none"> • Switch off instrument. Ensure adequate cooling, switch on again after prolonged cooling time.

7. Technical data

Properties	Values
General information	
Dimensions	mm
• Ring Light (OD x H, ID)	114 x 29, 66
• Transmitted Light Stage (OD x H)	84 x 16
• Spot Light Plus/Plus Double control unit (W x D x H)	~61 x 114 x 107
• Length flexible arm Spot light Plus / Plus Double incl. light head	500
Weight	kg
• Ring Light	~0,6
• Transmitted Light Stage	~0,7
• Spot Light Plus Single	~0,5
• Spot Light Plus Double	~0,7
Operating environment	°C +5 ... +40
Relative air humidity*	% up to 31°C ambient temperature: 85% up to 40°C ambient temperature: decreasing linearly to 75%
Air pressure*	hPa 800 ... 1060
Max. altitude	≤ 2000 m asl
Transport and storage	
• Temperature	°C -20 ... +70
• Rel. humidity	% 10 ... 95 (non condensing)
• Air pressure	hPa 500 ... 1200
Contamination level	2
IP-protection class	IP20 Following IEC 60529
Electrical information	
Operation voltage/ frequency	100–240 V ~ 50–60 Hz ±10 %
Nominal voltage	V
• Ring Light & Trans- mitted Light Stage	12 DC
• Spot Light	5 DC
Rated current	A max. 0.5
Power consumption	VA max. 25
Protection class	
• Power supply	II
• Lighting head	III
Overvoltage category	I
Lamp type	High Power LED
Average LED lifetime	
Setting level 5	50,000 (drop luminous flux to 70%)

Properties	Values	
Lighting specifications		
Max. illuminance	lx	
• Ring Light		90,000 (at 75 mm working distance)
• Ring Light Plus		140,000 (at 65 mm working distance)
Max. Luminance Transmitted Light Stage	cd/m ²	12,000
Total liminous flux at the single LED head of the spot lights (typ. value) level 5	lm	130
Farbtemperatur	K	approx. 5,600
Light exit angle		
Spot Light ($2\alpha_{\text{eff}}$)	°	approx. 60
Conformity		CE, UKCA (power supply) CE, UL, PSE
EMV emission class		B

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WEEE declaration

Your SCHOTT product was produced and developed with high quality materials and components. The symbol indicates that electrical and electronic devices must be separated from domestic waste and appropriately disposed of after useful life. SCHOTT AG Lighting and Imaging has arranged a waste management system for recycling. Please use this system for removal and help to protect the environment we live in.

Further information regarding our waste management system, please refer to www.schott.com/lightingimaging/recycle

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SCHOTT AG, Hattenbergstrasse 10,